

Rehabilitation bed type A-6-3S/T

- Four-segments, all segments are regulated electrically
- Wire remote controller
- Regulated backrest with auto-regress function
- Mattress base height adjusted electrically
- CPR/QR function of quick automatic adjustment to flat position of the whole bed
- Trendelenburg and anti-Trendelenburg positions
- Detachable foot and headboards rails
- Removable foot and headboards panels
- Different types of handrails available



Technical data :

| | |
|--|--------------------------|
| Length: total/mattress base | 2160 mm / 1900 – 2130 mm |
| Width with handrails/without handrails | 1 050 mm / 900 mm |
| Adjustment of the mattress base height | 400 ÷ 840 mm |
| Backrest segment angle adjustment | 0° ÷ 70° |
| Thigh rest segment angle adjustment | 0° ÷ 35° |
| Under thigh segment angle adjustment | - 5° ÷ 4° |
| Trendelenburg and anti-Trendelenburg pos. | 0 ÷ 25° ± 5 |
| Power requirements | ~ 230 V , 50Hz |
| Safety level | II B |



TRANSACTOR
Sp. z o.o.



ul. Sapiezynska 10
00-215 Warszawa
tel./fax: 22 635 59 02
e-mail: export-import@transactor.pl
www.transactor.pl