

Rehabilitation bed type B-1/3S

- Four segments , all segments regulated electrically
- Wire controller for all segments
- Electrical adjustment of bed's height
- Wooden-MDF bedstead with handrails
- Trendelenburg and anti-Trendelenburg position
- CPR/QR function – quick adjustment to flat position of the whole bed
- Individual lock for all wheels
- Detachable frame of the bed allows for easy transport in 4 pieces
- Easy assembly.

Technical data :

Length: total/mattress frame	2140 ± 5 mm / 2020 ± 5 mm
Width: total/mattress frame	1035 ± 5 mm / 890 ± 5 mm
Mattress frame height adjustment	400 ÷ 800 mm
Backrest segment angle adjustment	0° ÷ 75° ± 5°
Thigh section angle adjustment	0° ÷ 45° ± 5°
Under thigh segment angle	- 5° ÷ 4°
Trendelenburg and anti-Trendelenburg position	0 ÷ 25° ± 5°
Power requirement	~ 230 V , 50Hz
Safety level	II B



TRANSACTOR
Sp. z o.o.



ul. Sapiezynska 10
00-215 Warszawa
tel./fax: 22 635 59 02
e-mail: export-import@transactor.pl
www.transactor.pl